

Data Independent Acquisition (DIA)

Proteomic short Course 2025

Brett S. Phinney, Ph D.

UC Davis Genome Center

Before DIA

- Profiling was done by DDA
 - Spectral counting
 - LFQ
 - TMT (10 or 11 plex)

Example from a Paper we published in Plant Cell in 2019 using Spectral Counts and MS1 LFQ!

A	FET	Spectral counts					
		NC r1	NC r2	NC r3	ROG2 r1	ROG2 r2	ROG2 r3
Yeast/Mammals/Arabidopsis							
- /TRAPPC11/ROG2 (bait)	<0.0001	0	0	0	259	273	285
TRS20/TRAPPC2/At2g20930	<0.0001	0	0	0	12	6	14
TRS23/TRAPPC4/At5g02280	<0.0001	0	0	0	5	3	4
TRS31/TRAPPC5/At5g58030	<0.0001	0	0	0	7	9	12
TRS33/TRAPPC6/At3g05000	<0.0001	0	0	0	20	20	20
BET3/TRAPPC3/At5g54750	<0.0001	0	0	0	48	42	38
TRS85/TRAPPC8/At5g16280	<0.0001	0	0	0	92	85	81
TRS120/TRAPPC9/At5g11040	nd	nd	nd	nd	nd	nd	nd
TRS130/TRAPPC10/At5g54440	nd	nd	nd	nd	nd	nd	nd





Figure 7. AtTRAPPC11/ROG2 defines a novel Arabidopsis TRAPP complex.






(A) Spectral counts are presented for the six Arabidopsis TRAPP orthologs significantly enriched in YFP-ROG2 immunoprecipitates ("ROG2 r1-r3"), and for AtTRAPPC11/ROG2 itself (bait). r: biological replicate (1-3). FET: Fisher's exact test. NC: negative control. Note that Arabidopsis orthologs of TRS120 and TRS130, defining the TRAPP II complex, were not detected (nd) in our analysis. (B) Volcano plot illustrating log2 fold change (FC) (x axis) and statistical significance distribution (y axis) of the proteomic data set (LFQ intensities). Red squares indicate proteins significantly enriched in YFP-ROG2 immunoprecipitates. Analyses were performed in triplicates.

DDA vs DIA

DDA




DIA




In addition to the semi randomness of DDA. Peak picking is an issue with MS1 LFFQ

T:\Data\...\QEDDA\FL20190503_17_3-1_DDA


05/05/19 01:52:45



Typical? MS Spectra




Nice S/N





Typical? MS Spectra

FL20190503_17_3-1_DDA #12449 RT: 35.65 AV: 1 NL: 5.22E7
T: FTMS + p NSI Full ms [375.0000-1500.0000]




Nice S/N

FL20190503_17_3-1_DDA #12449 RT: 35.65 AV: 1 NL: 5.22E7
T: FTMS + p NSI Full ms [375.0000-1500.0000]




Lots of MS1 of peptides looks like this




FL20180720_37_HELA_UNIVERSAL #26236 RT: 63.98 AV: 1 NL: 5.00E7


T: FTMS + p NSI Full ms [375.0000-1500.0000]



FL20180720_37_HELA_UNIVERSAL #26236 RT: 63.98 AV: 1 NL: 1.52E7

T: FTMS + p NSI Full ms [375.0000-1500.0000]





Missing Values Comparison Peptide level

- DDA (maxquant 1% FDR protein and peptide)
- Fusion Lumos

Replicate	Missing Values
1	16.3%
2	15.15%
3	23.47%


Missing in all 3 = 5.53%

- In my hands Even with MBR on, it's difficult to get consistent peak picking across multiple samples


Real DIA example

Example MS1 Spectrum

FL20191105_51_HS2D2 #54209 RT: 42.04 AV: 1 NL: 2.18E8
T: FTMS + c NSI Full ms [390.0000-1010.0000]




Real Example MS1 Spectrum




1.5-2 sec total time

MS1 spectra zoomed in to around 750 m/z




T:\service\...\FL20191105_51_HS2D2

11/06/19 09:07:04




FL20191105_51_HS2D2 #54209 RT: 42.04 AV: 1 NI: 1.33E8


Real Example DIA spectrum 12 Da Window 750.5911 m/z



Next DIA spectrum 12 Da Window 762.5966 m/z




Remember LC-MS is three dimensional! You can use Retention time to check what fragment ions elute at the same time




T:\service\...\FL20191105_51_HS2D2

11/06/19 09:07:04




FL20191105_51_HS2D2 #54209 RT: 42.04 AV: 1 NI: 1.33E8


Example Peptide Centric detection



DIA Random peptide KSDGIYIINLK




DIA Random peptide KSDGIYIINLK



Why is the one on the right correct and the one on the left incorrect?

Check with DDA Data (DDA library)
Check with Deep learning (Prosit)



- DDA Library Spectra


bsa 100

13604.79785156

27.834156041.0000000029.31118774

U2002072309 4 (14.036) Sm (SG, 1x3.00); Sb (1,40.00); Cm (4:5)


2: TOF MSMS 733.28ES+
138



Many different ways to analyze DIA


- Spectrum Centric
 - Fragpipe
 - Reconstruct DDA spectrum from DIA Data by using RT
 - Match Spectra using traditional DDA methods
- Some say this scales better for large search spaces
- Peptide Centric
 - Spectronaut
 - DIA-NN
 - Tests every theoretical peptide in a proteome or subset of a proteome against the DIA Data to see if it's there.
 - Uses Deep learning to model theoretical peptide characterizes
 - RT
 - Fragment intensity
 - Fragment intensity ratios
 - It's all filtering and fragment / precursor correlation

Spectrum Centric






Makes a DDA spectra from ions correlated in time


Does PSM like traditional search engine against peptide sequence




Peptide Centric



Peptide Centric DIA




Ting 2015



Searle 2018


Precursor



— Precursor - - Precursor+1 isotope - - Precursor+2 isotope


☒ Show precursor isotopes

Fragments



— y3 — b4 — y4 — b5 — y5

DIA Random peptide KSDGIYIINLK



Chromatogram Libraries made with Deep learning!

nature > nature methods > articles > article



nature|methods

Article | Published: 27 May 2019

Prosit: proteome-wide prediction of peptide tandem mass spectra by deep learning


Siegfried Gessulat, Tobias Schmidt, Daniel Paul Zolg, Patroklos Samaras, Karsten Schnatbaum, Johannes Zerweck, Tobias Knaute, Julia Rechenberger, Bernard Delanghe, Andreas Huhmer, Ulf Reimer, Hans-Christian Ehrlich, Stephan Aiche, Bernhard Kuster✉ & Mathias Wilhelm✉

Using Deep Learning to Predict DDA libraries

One issue with peptide DIA is that it seems to work a lot better with a DDA Library

- You need to generate one yourself
- Or use a public one (Pan Human)
- But now, with **Prosit**, you can use deep learning to predict one without generating one using your mass spec!
- Introduces at ASMS this year

Fig. 1: Accurate retention time and fragment ion intensity prediction by deep learning.



Article | Published: 27 May 2019

Prosit: proteome-wide prediction of peptide tandem mass spectra by deep learning

Siegfried Gessulat, Tobias Schmidt, Daniel Paul Zolg, Patroklos Samaras, Karsten Schnatbaum, Johannes Zerweck, Tobias Knaute, Julia Rechenberger, Bernard Delanghe, Andreas Huhmer, Ulf Reimer, Hans-Christian Ehrlich, Stephan Aiche, Bernhard Kuster & Mathias Wilhelm

Nature Methods **16**, 509–518 (2019) | [Download Citation](#)



Cold Spring Harbor Laboratory

bioRxiv
THE PREPRINT SERVER FOR BIOLOGY

HOME | A

Search

New Results

[Comment on this paper](#)

Generating high-quality libraries for DIA-MS with empirically-corrected peptide predictions

Brian C. Searle, Kristian E. Swearingen, Christopher A. Barnes, Tobias Schmidt, Siegfried Gessulat, Bernhard Kuster, Mathias Wilhelm

doi: <https://doi.org/10.1101/682245>



This article is a preprint and has not been peer-reviewed [what does this mean?].

Prosit used millions of synthetic peptides to build deep learning models


- Predicts what y and b ion you should see
- Predicts the intensity of fragment ions
- Predicts retention time

Pros: It works!


Cons: does not tell us anything about the underlying chemistry



These intensities can now be predicted accurately




We Can check if these Ratios are right!



We can also use deep learning to check if the RT is right!

DIA Random peptide KSDGIYIINLK



Why is the one on the right matches the predicted ratios the one on the left does not